



APOLLO

IMPLANT COMPONENTS



SCREW RETAINED CASES

APOLLO IMPLANT COMPONENTS created a special offer for multi-implant cases: Structures, Supra-Structures, Implant Bridges and Implant Bars based on multi-unit abutments. Thanks to this solution each case has a perfect passivity and tightness. The offer is available for most implant systems.

APOLLO SCREW RETAINED CASES																		
IMPLANT SYSTEM compatible with:	PLATFORM	CUSTOM MULTIUNIT	MULTIUNIT LEVEL (MU)		IMPLANT LEVEL (MP)		ANGULATED SCREW CHANNEL up to 28° *		DIGITAL WORKFLOW				PROSTHETICS SCREW					
			MATERIAL		MATERIAL		MATERIAL		use APOLLO SCANBODY marked as:		choose dedicated CAD library component (MP/MU):		REF. numbers			Torque [Ncm]	Standard Screwdriver	
									multiunit level (MU)	implant level (MP)	for EXOCAD		for 3SHAPE		APOLLO screw			DYNAMIC screw**
			Ti	CoCr	Ti	CoCr	Ti	CoCr	SCANBODY MULTIUNIT	SCANBODY GEN II								APOLLO / Dynamic
AB	3.5	x			x	x	x	x		AAB 35	AAB 3.5 MP Non Engaging GEN II Scanbody	AAB 3.5 MP Non Engaging	005-038-01-03	41.317.071.01-2	D4A-S	25/25 Ncm	1.25 / HEX	
Adin Touareg Close Fit	Narrow	x			x	x				ADCLF N	ADCLF N MP Non Engaging GEN II Scanbody	ADCLF N MP Non Engaging	005-002-01	-	NP0001	35/- Ncm	1.27 / HEX	
	Regular	x			x	x	x			ADCLF R	ADCLF R MP Non Engaging GEN II Scanbody	ADCLF R MP Non Engaging	005-002-02	41.316.073.01-2	RP0001	35/20 Ncm		
	Wide				x	x	x	x		ADCLF W	ADCLF W MP Non Engaging GEN II Scanbody	ADCLF W MP Non Engaging	005-002-03	41.320.075.01-2	WP0001	35/25 Ncm		
Adin Touareg S & OS	3.5	x			x	x	x	x		ADSOS	ADSOS MP Non Engaging GEN II Scanbody	ADSOS MP Non Engaging	005-038-01-03	-	RS3402	25/- Ncm	1.27 / HEX	
AlphaBio SPI	3.75 - 5.3	x			x	x	x	x		ALPHB SPI	ALPHB SPI MP Non Engaging GEN II Scanbody	ALPHB SPI MP Non Engaging	005-038-01-03	41.317.071.01-2	5122	30/25 Ncm	1.25 / HEX	
	Multi-unit «4,7		x	x					ALPHB MU	ALPHB MU Non Engaging	ALPHB MU Non Engaging	005-004-M47	-	6092	25/- Ncm	1.25 / HEX		
AlphaBio Nice	3.2				x	x	x		ALPHB MU	ALPHB Nce	ALPHB Nce MP Non Engaging GEN II Scanbody	ALPHB Nce MP Non Engaging	005-005-01	-	7345	20/- Ncm	1.25 / HEX	
	Multi-unit «4,7		x	x						ALPHB MU Non Engaging	ALPHB MU Non Engaging	005-004-M47	-	6092		1.25 / HEX		
Ankylos	C	x			x	x			ANK CX	ANK C MP Non Engaging GEN II Scanbody	ANK C MP Non Engaging	005-006-01-02 (integrated)	-	-	15/- Ncm	1.00 / HEX		
	CX	x			x	x				ANK CX MP Non Engaging GEN II Scanbody	ANK CX MP Non Engaging	-	-	15/- Ncm				
	Multi-unit «4,2		x	x					ANK MU 42	ANK MU 4.2 Non Engaging	ANK MU 4.2 Non Engaging	005-006-M42	-	31056022	15/- Ncm			
	Multi-unit «5,5		x	x					ANK MU 55	ANK MU 5.5 Non Engaging	ANK MU 5.5 Non Engaging	005-006-M55	-	31056021	15/- Ncm			
Megagen AnyRidge	3.5-8.0	x			x	x	x	x		MGANY	MGANY MP Non Engaging GEN II Scanbody	MGANY MP Non Engaging	005-007-01	41.318.075.01-2	AANM	30/25 Ncm	1.20 / HEX	
Astra OsseoSpeed EV	3.0				x	x				AEV 30	AEV 3.0 MP Non Engaging GEN II Scanbody	AEV 3.0 MP Non Engaging	005-008-01	-	25203	25/- Ncm	1.25 / HEX	
	3.6				x	x	x			AEV 36	AEV 3.6 MP Non Engaging GEN II Scanbody	AEV 3.6 MP Non Engaging	005-008-02	41.316.072.01-2	25204	25/20 Ncm		
	4.2				x	x	x	x		AEV 42	AEV 4.2 MP Non Engaging GEN II Scanbody	AEV 4.2 MP Non Engaging	005-008-03	41.318.074.01-2	25205	25/25 Ncm		
	4.8				x	x	x	x		AEV 48	AEV 4.8 MP Non Engaging GEN II Scanbody	AEV 4.8 MP Non Engaging	005-008-04-05	41.320.074.01-2	25206			
	5.4				x	x	x	x		AEV 54	AEV 5.4 MP Non Engaging GEN II Scanbody	AEV 5.4 MP Non Engaging		41.320.074.01-2	25207			
	Multi-unit «3,7		x	x			x	x	AEV MU	AEV MU Non Engaging	AEV MU Non Engaging	005-008-M37	41.318.045.01-2	25481	25/25 Ncm			
Astra OsseoSpeed TX	Yellow 3.0	x			x	x			ATX MU	ATX 30	ATX 3.0 MP Non Engaging GEN II Scanbody	ATX 3.0 MP Non Engaging	005-009-01	-	24289	15/- Ncm	1.25 / HEX	
	Aqua 3.5/4.0	x			x	x	x			ATX 35/40	ATX 3.5/4.0 MP Non Engaging GEN II Scanbody	ATX 3.5/4.0 MP Non Engaging	005-009-02	41.316.076.01-2	24449	20/20 Ncm		
	Lilac 4.5/5.0				x	x	x	x		ATX 45/50	ATX 4.5/5.0 MP Non Engaging GEN II Scanbody	ATX 4.5/5.0 MP Non Engaging	005-009-03	41.320.090.01-2	24209	25/25 Ncm		
	Multi-unit «3,4		x	x						ATX MU Non Engaging	ATX MU Non Engaging	005-009-M34		22435	15/- Ncm			
Biomet 3i Certain	3.4	x			x	x	x		B3IC 41/50	B3IC 34	B3IC 3.4 MP Non Engaging GEN II Scanbody	B3IC 3.4 MP Non Engaging	005-011-01-03	41.316.084.01-2	IUNIHG	20/20 Ncm	1.25 / HEX	
	4.1	x			x	x	x			B3IC 4.1 MP Non Engaging GEN II Scanbody	B3IC 4.1 MP Non Engaging							
	5.0	x			x	x	x			B3IC 5.0 MP Non Engaging GEN II Scanbody	B3IC 5.0 MP Non Engaging							
Biohorizons	3.0	x			x	x			BHZON 57	BHZON 30	BHZON 3.0 MP Non Engaging GEN II Scanbody	BHZON 3.0 MP Non Engaging	005-012-01-04	41.317.071.01-2	PXAS	30/- Ncm	1.25 / HEX	
	3.5	x			x	x	x	x		BHZON 35	BHZON 3.5 MP Non Engaging GEN II Scanbody	BHZON 3.5 MP Non Engaging				30/25 Ncm		
	4.5				x	x	x	x		BHZON 45	BHZON 4.5 MP Non Engaging GEN II Scanbody	BHZON 4.5 MP Non Engaging						
	5.7				x	x	x	x		BHZON 57	BHZON 5.7 MP Non Engaging GEN II Scanbody	BHZON 5.7 MP Non Engaging						
	Multi-unit «4,8		x	x					BHZON MU	BHZON MU Non Engaging	BHZON MU Non Engaging	005-012-M48		PXMUPSR / PXMUPSS	15/- Ncm			
Bredent Sky	3.5				X	X	X	X	BRESK	BRESK 3.5 MP Non Engaging GEN II Scanbody	BRESK 3.5 MP Non Engaging	005-013-01-04	41.318.083.01-2	SKY-PS22	25/25 Ncm	TORX		
	4,0				X	X	X	X		BRESK 4.0 MP Non Engaging GEN II Scanbody	BRESK 4.0 MP Non Engaging							
	4,5				X	X	X	X		BRESK 4.5 MP Non Engaging GEN II Scanbody	BRESK 4.5 MP Non Engaging							
	5,5				X	X	X	X		BRESK 5.5 MP Non Engaging GEN II Scanbody	BRESK 5.5 MP Non Engaging							
CWM INNO-SUBMERGED	Narrow	x			x	x	x		CWIS R	CWIS N MP Non Engaging GEN II Scanbody	CWIS N MP Non Engaging	005-031-01	41.316.094.01-2	SSHR100N	20/20 Ncm	1.20 / HEX		
	Regular	x			x	x	x	x		CWIS R MP Non Engaging GEN II Scanbody	CWIS R MP Non Engaging	005-025-01-02	41.320.079.01-2	2SSHR200	30/25 Ncm			
Dentaurum Tiologic	S 3.3				x	x	x		TIOL S	TIOL S MP Non Engaging GEN II Scanbody	TIOL S MP Non Engaging	005-017-01-03	41.316.081.01-2	386-310-00	30/20 Ncm	1.30 / HEX		
	M (3.7-4.2)				x	x	x			TIOL M MP Non Engaging GEN II Scanbody	TIOL M MP Non Engaging							
	L (4.8 -5.5)				x	x	x			TIOL L MP Non Engaging GEN II Scanbody	TIOL L MP Non Engaging							
Dentium SuperLine	3.6 - 7.0	x			x	x	x	x	DENS	DENS MP Non Engaging GEN II Scanbody	DENS MP Non Engaging	005-018-01	41.320.079.01-2	ASC2045	30/25 Ncm	1.28 / HEX		
	Multi-unit «4,3		x	x				DENS MU 43		DENS MU 4.3 Non Engaging	DENS MU 4.3 Non Engaging	005-018-M43	-	GSRS16T	20/- Ncm			
	Multi-unit «4,5		x	x					DENS MU 45/55	DENS MU 4.5 Non Engaging	DENS MU 4.5 Non Engaging	005-018-M45-M55	-	SRS18T	30/- Ncm			
	Multi-unit «5,5		x	x						DENS MU 5.5 Non Engaging	DENS MU 5.5 Non Engaging		-					
GlobalD In-Kone	3.5-5.0	x			x	x			GDINK	GDINK MP Non Engaging GEN II Scanbody	GDINK MP Non Engaging	005-019-01	-	DVPCI	15/- Ncm	1.20 / HEX		
	Multi-unit «4,3		x	x						Original GlobalD Scanbody/Library		005-019-M43	-	DVPIINK	15/- Ncm			
ICX Templant	3.75 - 4.8	x			x	x			ICXTMP MU	ICXTmp MP Non Engaging GEN II Scanbody	ICXTmp MP Non Engaging	005-020-01	-	C-007-000001-01	30/- Ncm	1.40 / HEX		
	Multi-unit «4,8		x	x						ICXTmp MU Non Engaging	ICXTmp MU Non Engaging	005-020-M48	-	C-020-000011-01	27/- Ncm	1.20 / HEX		
Implant Direct Legacy	3.5	x			x	x	x	x	IMPLCY 35	IMPLCY 3.5 MP Non Engaging GEN II Scanbody	IMPLCY 3.5 MP Non Engaging	005-038-01-03	41.317.071.01-2	1000-01	25/25 Ncm	1.25 / HEX		
MIS	Multi-unit «4,8		x	x			x		Original MIS Scanbody/Library		005-022-M48	41.316.044.01-2	MU-S0220	30/20 Ncm	1.27 / HEX			
Medentika MedentIBASE	Multi-unit «3,7		x	x					Original MEDENTIKA Scanbody/Library			005-041-M37	-	B7600	15/- Ncm			
NeoBiotech	3.5	x			x	x	x	x	NEOB MU	NEOB 3.5 MP Non Engaging GEN II Scanbody	NEOB 3.5 MP Non Engaging	005-025-01-02	41.320.079.01-2	ISC520	25/25 Ncm	1.20 / HEX		
	4.0 - 7.0	x			x	x	x	x		NEOB MP Non Engaging GEN II Scanbody	NEOB MP Non Engaging							
	Multi-unit «4,8		x	x						NEOB MU Non Engaging	NEOB MU Non Engaging	005-025-M48	-	MUAS14	15/- Ncm			
Neodent CM	3.5 - 5.0	x			x	x			NDENT CM	NDENT CM MP Non Engaging GEN II Scanbody	NDENT CM MP Non Engaging	005-006-01-02 (int)	-	116038	15/- Ncm	1.00 / HEX		
Neodent GM	3.5 - 5.0	x			x	x	x		NDENT CM	NDENT CM MP Non Engaging GEN II Scanbody	NDENT CM MP Non Engaging	005-011-01-03	41.316.084.01-2	116286	20/20 Ncm	1.25 / HEX		
Nobel Biocare Active CC	3.0	x			x	x			NACT NP	NACT 30	NACT 3.0 MP Non Engaging GEN II Scanbody	NACT 3.0 MP Non Engaging	005-028-01	-	37890	15/- Ncm	UNIGRIP	
	NP 3.5	x			x	x	x			NACT NP MP Non Engaging GEN II Scanbody	NACT NP MP Non Engaging	005-028-02	41.316.073.01-2	37891	35/20 Ncm			
	RP 3.4/5.0	x			x	x	x	x		NACT RP MP Non Engaging GEN II Scanbody	NACT RP MP Non Engaging	005-028-03	41.320.075.01-2	37892	35/25 Ncm			
	Multi-unit «4,8		x	x			x		NOB MU	NOB MU Non Engaging	NOB MU Non Engaging	005-028-M48	41.314.039.01-2	29285	15/20 Ncm			
Nobel Biocare Replace	NP 3.5				x	x	x	x	NRPL 55	NRPL 35	NRPL 3.5 MP Non Engaging GEN II Scanbody	NRPL 3.5 MP Non Engaging	005-030-01	41.318.075.01-2	36818	35/25 Ncm		
	RP 4.3/5.0				x	x	x	x		NRPL 43	NRPL 4.3 MP Non Engaging GEN II Scanbody	NRPL 4.3 MP Non Engaging	005-030-02	41.320.090.01-2	29475			
	WP 5.5				x	x	x	x		NRPL 55	NRPL 4.3 MP Non Engaging GEN II Scanbody	NRPL 4.3 MP Non Engaging	005-030-03	-				
	Multi-unit «4,8		x	x			x		NOB MU	NOB MU Non Engaging	NOB MU Non Engaging	005-028-M48	41.314.039.01-2	29285	15/20 Ncm			
Osstem TS	Mini 3.5	x			x	x	x		OSSTS M	OSSTS M MP Non Engaging GEN II Scanbody	OSSTS M MP Non Engaging	005-031-01	41.316.094.01-2	GSAB5M	20/20 Ncm	1.20 / HEX		
	Regular (4.0 -7.0)	x			x	x	x	x		OSSTS R	OSSTS R MP Non Engaging GEN II Scanbody	OSSTS R MP Non Engaging	005-031-02	41.320.079.01-2	GSABSS		30/25 Ncm	
	Multi-unit «4,8		x	x					OSSTS MU	OSSTS MU Non Engaging	OSSTS MU Non Engaging	005-031-M48	-	MTS200	20/- Ncm			
	Convertible 4.0		x	x					OSSTS CNV	OSSTS CNV 4.0 Non Engaging	OSSTS CNV 4.0 Non Engaging	005-031-M40C	-	GSFSM	20/- Ncm			
Southern Implants External Hex	3.0				x	x			SIEXH 375-50	SIEXH 30	SIEXH 3.0 MP Non Engaging GEN II Scanbody	SIEXH 3.0 MP Non Engaging	005-033-01		TS-P-16	25-32/- Ncm	1.22 / HEX	
	3.25				x	x				SIEXH 325	SIEXH 3.25 MP Non Engaging GEN II Scanbody	SIEXH 3.25 MP Non Engaging		-				
	3.75/4,0				x	x				SIEXH 375/4,0 MP Non Engaging GEN II Scanbody	SIEXH 3.75/4.0 MP Non Engaging	005-033-02-04	-	TSH3				
	5.0				x	x				SIEXH 5.0 MP Non Engaging GEN II Scanbody	SIEXH 5.0 MP Non Engaging		-					
Straumann Bone Level	SC 2.9	x			x	x			STRBL 41	STRBL 29	STRBL 2.9 MP Non Engaging GEN II Scanbody	STRBL 2.9 MP Non Engaging	005-034-01	-	025.0031	35/- Ncm		
	N.C 3.3	x			x	x	x			STRBL 33	STRBL 3.3 MP Non Engaging GEN II Scanbody	STRBL 3.3 MP Non Engaging	005-034-02-03	41.316.078.01-2	025.2900	35/20 Ncm		
	R.C 4.1	x			x	x	x		STRBL 41	STRBL 4.1 MP Non Engaging GEN II Scanbody	STRBL 4.1 MP Non Engaging	025.4900						
	Multi-unit «3,5		x	x					STR MU 35	STR MU 3.5 Non Engaging	STR MU 3.5 Non Engaging	005-034-M35	-	023.4763	35/- Ncm			
	Multi-unit «4,6		x	x			x		STR MU 46	STR MU 4.6 Non Engaging	STR MU 4.6 Non Engaging	005-034-M46	-	023.4763				
Straumann Tissue Level	NNC 3.5				x	x			STRTL 48	STRTL NNC	STR							



APOLLO

IMPLANT COMPONENTS



APOLLO SCREW RETAINED CASES																			
IMPLANT SYSTEM compatible with:	PLATFORM	CUSTOM MULTIUNIT	MULTIUNIT LEVEL (MU)		IMPLANT LEVEL (MP)		ANGULATED SCREW CHANNEL up to 28° *		DIGITAL WORKFLOW				PROSTHETICS SCREW						
									use APOLLO SCANBODY marked as:		choose dedicated CAD library component (MP/MU):		REF. numbers			Torque [Ncm]	Standard Screwdriver		
			MATERIAL		MATERIAL		MATERIAL		multiunit level (MU)	implant level (MP)	for EXOCAD		for 3SHAPE		APOLLO screw			DYNAMIC screw**	Implant manufacturer screw
			Ti	CoCr	Ti	CoCr	Ti	CoCr	SCANBODY MULTIUNIT	SCANBODY GEN II									
Straumann BLX	RB	x			x	x			STRX RB	STRX RB MP Non Engaging GEN II Scanbody		STRX RB MP Non Engaging	005-036-01-02	-	065.0036	35/- Ncm			
	WB				x	x				STRX WB MP Non Engaging GEN II Scanbody		STRX WB MP Non Engaging		-					
	Multi-unit e3.5		x	x				STR MU 35	STR MU 3.5 Non Engaging		STR MU 3.5 Non Engaging	005-034-M35	-	023.4763	35/- Ncm				
	Multi-unit e4.6		x	x			x	STR MU 46	STR MU 4.6 Non Engaging		STR MU 4.6 Non Engaging	005-034-M46	-	023.4763					
Thommen	3.5				x				THO 35	THO 3.5 MP Non Engaging GEN II Scanbody		THO 3.5 MP Non Engaging	005-037-01	-	403500	15/- Ncm	4-LOBE		
	4.0				x		x		THO 40	THO 4.0 MP Non Engaging GEN II Scanbody		THO 4.0 MP Non Engaging	005-037-02	41.316.076.01-2	403501	25/20 Ncm			
	4.5				x		x		THO 45	THO 4.5 MP Non Engaging GEN II Scanbody		THO 4.5 MP Non Engaging							
	5.0				x		x			THO 5.0 MP Non Engaging GEN II Scanbody		THO 5.0 MP Non Engaging							
	Multi-unit e4.8		x	x				THO MU	THO MU Non Engaging		THO MU Non Engaging	005-037-M48	-	4.03.507	15/- Ncm				
Zimmer Tapered Screw Vent	3.5	x			x	x	x	x	ZIMTS 35	ZIMTS 3.5 MP Non Engaging GEN II Scanbody		ZIMTS 3.5 MP Non Engaging	005-038-01-03	41.317.071.01-2	MHLAS	30/25 Ncm	1.25 / HEX		
	4.5				x	x	x	x	ZIMTS 45	ZIMTS 4.5 MP Non Engaging GEN II Scanbody		ZIMTS 4.5 MP Non Engaging							
	5.7				x	x	x	x	ZIMTS 57	ZIMTS 5.7 MP Non Engaging GEN II Scanbody		ZIMTS 5.7 MP Non Engaging							
Zimmer Eztec	3.1	x			x	x			ZIMEZ 31	ZIMEZ 3.1 MP Non Engaging GEN II Scanbody		ZIMEZ 3.1 MP Non Engaging	005-039-01	-	CKAS	30/- Ncm			
XIVE	3.0	x			x	x	x		XVE 30	XVE 3.0 MP Non Engaging GEN II Scanbody		XVE 3.0 MP Non Engaging	005-040-01	-	32464304	24/- Ncm		1.22 / HEX	
	3.4	x			x	x	x		XVE 34	XVE 3.4 MP Non Engaging GEN II Scanbody		XVE 3.4 MP Non Engaging	005-040-02	41.316.081.01-2	32454305	24/20 Ncm			
	3.8	x			x	x	x		XVE 38	XVE 3.8 MP Non Engaging GEN II Scanbody		XVE 3.8 MP Non Engaging							
	Multi-unit e4.8		x	x				XVE MU 48		XVE MU 4.8 Non Engaging		XVE MU 4.8 Non Engaging	005-040-M48	-	25-4300	24/- Ncm			
	Multi-unit e5.2		x	x						XVE MU 5.2 Non Engaging		XVE MU 5.2 Non Engaging	005-040-M52	-	25-4300	24/- Ncm			

All trademarks and company names are the property of their respective owners.

* Screw channel angulation option excludes choosing retention.

** Screw and screwdriver produced by Dynamic Abutment® Solutions.

x - available

Document update at 15.04.2020.

ORDER VIA APOLLO ORDER PORTAL
www.apollocomponents.eu

3shape

APOLLO LIBRARY
officially available on the 3Shape Download Center library list.

exocad

APOLLO LIBRARY
officially available on the exocad DentalCAD library list.

ANGULATED SCREW-HOLE SOLUTIONS

The angulated screw-hole provides the option to place the screw access hole anywhere between 0° to 28° in a 360° radius.

AUTHORIZED MILLING CENTER
DYNAMIC ABUTMENT® SOLUTIONS



HELP DESK



ORDER VIA APOLLO ORDER PORTAL
www.apollocomponents.eu



+48 42 225 29 29



support@apollocomponents.eu



POLAND | 95-200 Pabianice | Konopna 16